

## EasyPfu DNA Polymerase

**Cat. No.** AP211

**Concentration** 2.5 units/ $\mu$ l

**Storage** -20°C for two years

### Description

*EasyPfu* DNA Polymerase is an engineered version of *pfu* DNA Polymerase with enhanced yield and higher fidelity. *EasyPfu* DNAPolymerase possesses a proofreading 3'-5' exonuclease activity.

### Highlights

- *EasyPfu* DNA Polymerase offers 18-fold fidelity as compared to *EasyTaq*® DNA Polymerase.
- Extension rate is about 0.5 kb/min.
- PCR products can be directly cloned into *pEASY*®-Blunt vectors.
- Amplification of genomic DNA fragment up to 6 kb.
- Amplification of plasmid DNA fragment up to 10 kb.

### Applications

- High fidelity PCR
- Blunt-end cloning
- Site-directed mutagenesis

### Unit Definition

One unit of *EasyPfu* DNA Polymerase incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

### Quality Control

*EasyPfu* DNA Polymerase has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity, >99% homogeneous measured by SDS-PAGE. Each batch of *EasyPfu* DNA Polymerase has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.

### Storage Buffer

50 mM Tris-HCl (pH 8.0), 50 mM KCl, 1 mM DTT, 0.1 mM EDTA, 50% (v/v) glycerol stabilizers

### 10×*EasyPfu* Buffer with 20 mM MgSO<sub>4</sub>

200 mM Tris-HCl (pH 8.8), 100 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 100 mM KCl, 20 mM MgSO<sub>4</sub>, others

### Kit Contents

Component	AP211-01/11	AP211-02/12	AP211-03/13
<i>EasyPfu</i> DNA Polymerase	250 U×1	500 U×1	500 U×6
10× <i>EasyPfu</i> Buffer	1.2 ml ×1	1.2 ml ×1	1.2 ml ×6
2.5 mM dNTPs	- / 500 $\mu$ l ×1	- / 1 ml ×1	- / 1 ml ×6
6×DNA Loading Buffer	500 $\mu$ l×1	1 ml ×1	1 ml ×2
50 mM MgSO <sub>4</sub>	200 $\mu$ l×1	400 $\mu$ l×1	1 ml ×1

### Reaction Components

Component	Volume	Final Concentration
Template	Variable	as required
Forward Primer (10 µM)	1 µl	0.2 µM
Reverse Primer (10 µM)	1 µl	0.2 µM
10× <i>EasyPfu</i> Buffer	4 µl	1×
2.5 mM dNTPs	5 µl	0.25 mM
<i>EasyPfu</i> DNA Polymerase	1 µl	2.5 units
ddH <sub>2</sub> O	Variable	-
Total volume	50 µl	-

### Thermal cycling conditions

94°C	2-5 min	} 30-35 cycles
94°C	30 sec	
50-60°C	30 sec	
72°C	0.5 kb/min	
72°C	5-10 min	

### Notes

- For GC-rich templates, the recommended denaturation temperature is 98°C.
- To ensure high fidelity, we recommended using high quality dNTPs. dNTPs containing dUTP cannot be used.
- Since it is not hot-start, we recommended to add enzyme last during PCR.

**FOR RESEARCH USE ONLY**