

## FavorPrep™ Total RNA Plus Mini Kit

- For isolation total RNA from animal cells and tissues.
- Efficient removal of genomic DNA by using gDNA Removal Column, without the need for DNase I digestion.

### Kit Contents:

For Research Use Only

Cat. No:	FATRK-P-004 (4 preps)	FATRK-P-050 (50 preps)	FATRK-P-100 (100 preps)
Lysis Buffer RXB	1.6 ml	20 ml	40 ml
Wash Buffer (concentrate) *	1.5 ml	15 ml	35 ml
RNase-free Water	0.5 ml	6 ml	6 ml
gDNA Removal Column (green)	4 pcs	50 pcs	100 pcs
RNA Mini Column	4 pcs	50 pcs	100 pcs
Collection Tube	12 pcs	150 pcs	300 pcs
Elution Tube	4 pcs	50 pcs	100 pcs
User Manual	1	1	1
Preparation of Wash Buffer by adding ethanol (96 ~ 100%)			
* Ethanol volume for Wash Buffer	6 ml	60 ml	140 ml

### Specification:

Principle: spin column (silica membrane)


Sample size : animal cells up to  $10^7$   
Tissue up to 30 mg

Elution volume: 30 ~ 50  $\mu$ l

### Important Notes:

1. Make sure the starting sample amount is under the limit.
2. Make sure everything is RNase-free when handling RNA.
3. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
4. Add required volume of RNase-free ethanol (96~100%) to Wash Buffer when first use.

### Brief Procedure:



## Protocol: Isolation of Total RNA from Animal Cells

Please Read Important Notes Before Starting Following Steps.

Additional requirement:  $\beta$ -Mercaptoethanol  
70% RNase-free ethanol  
rotor-stator homogenizer or 20-G needle syringe

1. Collect up to  $10^7$  cells by centrifuge at  $300 \times g$  for 5 min at  $4^\circ\text{C}$ . Remove all the supernatant.  
Add  $350 \mu\text{l}$  of Lysis Buffer RXB and  $3.5 \mu\text{l}$  of  $\beta$ -Mercaptoethanol to the cell pellet. Vortex vigorously for 1 min to resuspend the cells completely.  
**-- Note! Do not overload, too much sample will make cell lysis incompletely and lead to lower RNA yield and purity.**
2. Homogenize the sample by using a rotor-stator homogenizer or by passing the sample lysate through a 20-G needle syringe 10 times.  
**-- Important step: In order to release more RNA from samples, it is recommended to homogenize the sample by using suitable homogenize equipment, for example, with a rotor-stator homogenizer.**
3. Incubate at room temperature for 5 min.
4. Place a gDNA Removal Column to a Collection Tube and transfer the sample lysate to the gDNA Removal Column.
5. Centrifuge at full speed ( $\sim 18,000 \times g$ ) for 1 min. After centrifugation, do not discard the flow-through inside the Collection tube.
6. Transfer the supernatant of flow-through from step 5 to a 1.5 ml tube (not provided). Measure the volume of the supernatant.
7. Add 1 volume of 70% ethanol and mix well by plus-vortexing.  
Example: add  $330 \mu\text{l}$  of 70 % ethanol to  $330 \mu\text{l}$  of supernatant from step 6.
8. Place a RNA Mini Column in a Collection Tube and transfer the sample mixture to the RNA Mini Column.
9. Centrifuge at full speed ( $\sim 18,000 \times g$ ) for 1 min. Discard the flow-through and place the RNA Mini Column back to Collection Tube.
10. Add  $500 \mu\text{l}$  of Wash Buffer to the RNA Mini Column. Centrifuge at full speed ( $\sim 18,000 \times g$ ) for 1 min. Discard the flow-through and return the RNA Mini Column back to the Collection Tube.  
**-- Note: Make sure that ethanol has been added to Wash Buffer when first use.**
11. Repeat step 10 for one more washing.
12. Centrifuge the RNA Mini Column at full speed ( $\sim 18,000 \times g$ ) for an additional 3 min to dry the RNA Mini Column.  
**-- Important Step! This step will avoid the residual liquid to inhibit subsequent enzymatic reaction.**
13. Place the RNA Mini Column to a Elution Tube (provided, 1.5 ml microcentrifuge tube).
14. Add  $30 \sim 50 \mu\text{l}$  of RNase-free ddH<sub>2</sub>O to the membrane center of the RNA Mini Column. Stand the RNA Mini Column at room temperature for 1 min.  
**-- Important Step! For effective elution, make sure that RNase-free ddH<sub>2</sub>O is dispensed on the membrane center and is absorbed completely.**  
**-- Important : Do not elute the RNA using RNase-free water less than suggested volume. It will lower the RNA yield.**
15. Centrifuge the RNA Mini Column at full speed ( $\sim 18,000 \times g$ ) for 1 min to elute RNA. Store RNA at  $-80^\circ\text{C}$ .

## Protocol: Isolation of Total RNA from Animal Tissues

Please Read Important Notes Before Starting Following Steps.

Additional equipment: liquid nitrogen & mortar  
a rotor-stator homogenizer or a 20-G needle syringe  
 $\beta$ -Mercaptoethanol  
70% RNase-free ethanol

1. Weight up to 30 mg of tissue sample. Grind the sample in liquid nitrogen to a fine powder with a mortar and transfer the powder to a new microcentrifuge tube (not provided). Add  $350 \mu\text{l}$  of Lysis Buffer RX and  $3.5 \mu\text{l}$  of  $\beta$ -Mercaptoethanol.
2. Homogenize the sample by using a rotor-stator homogenizer or by passing the sample lysate through a 20-G needle syringe 10 times.  
**-- Note! Avoid thawing the sample during weighing and grinding.**  
**-- Note! Do not overload, too much sample will make cell lysis incompletely and lead to lower RNA yield and purity.**
3. Follow the Animal Cells Protocol starting from step 3.