



www.civicbio.com
info@civicbio.com orders@civicbio.com 1-877-661-4711

Livraison Gratuite
commandes > 100\$
25\$ si < 100\$

Récompenses Cashback

Produits Illimités

Promos Printemps 2018 !

Fin : 30 juin 2018

Elabscience®


2 pour 1 **ANTICORPS PRIMAIRES**
20, 60, 120, 200 ul



Sur tous les

TRANS®

2 pour 1 EvaGreen qPCR,
DNA ladders,
Competent cells
Cloning Kits
etc

NOS MEILLEURS VENDEURS! !
15 - 30%




 **Cashback x2 Promotion** 

DNA/RNA Purification & Extraction Kits

Micropipette ou Rack Magnétique **GRATUIT!**

PRODUITS & SERVICES
CIVIC BIOSCIENCE
 LE CIVISME AU CŒUR DE VOTRE RECHERCHE



www.civicbio.com
 info@civicbio.com orders@civicbio.com 1-877-661-4711

Livraison Gratuite

commandes > 100\$
 25\$ si < 100\$

Récompenses Cashback

Produits Illimités

Promos Printemps 2018

Fin : 12 juin 2018

Elabscience®
2 pour 1

Sur tous les
**ANTICORPS
 PRIMAIRES**

20, 60, 120, 200 ul



Qualité Garantie

Réactivité et Applications validées
 (WB, IHC, IF, IP, FC, etc)




Beaucoup de Choix

> 10 000 anticorps mono-
 et polyclonaux disponibles




Formats variés

Volumes de 20, 60, 120 et 200 ul.



Western Blot analysis of different cell lysates using PKC gamma Monoclonal Antibody E-AB-27770 at 1:5000 dilution.



Immunohistochemistry of paraffin-embedded Human colon cancer using PPP1R1B Polyclonal Antibody - E-AB-14928.

Elabscience®

~ 5000 Kits

ELISA / CLIA

-15% sur TOUT


Prix de volume aussi disponible!

Réactivité

Universal	Chicken
Human	Guinea Pig
Mouse	Monkey
Rat	Pig
Rabbit	Sheep

Populaires

Folic Acid E-EL-0009
 Vitamin D E-EL-0012
 Mouse TNF-α E-EL-M0049
 Human IL-6 E-EL-H0102
 Human C3a E-EL-H0818
 Human MMP-10 E-EL-H0132
 Rat Insulin E-EL-R2466



Livraison Gratuite

commandes > 100\$
 25\$ si < 100\$

Récompenses Cashback

Produits Illimités

Promos Printemps 2018

Fin : 30 juin 2018 ou 31 déc. 2018 pour le qPCR!

2 pour 1

TransStart TIP
EvaGreen
 qPCR
 SuperMix

qPCR mix comparison

	Gene	qPCR mix	Ct
StepOne Plus	18S	TransStart® TIP Green	18.52
		TransStart® TOP Green	18.07
		Power SYBR® Green	22.90
	β-actin	TransStart® TIP Green	12.39
		TransStart® TOP Green	13.91
		Power SYBR® Green	16.98
	GAPDH	TransStart® TIP Green	23.12
		TransStart® TOP Green	23.73
		Power SYBR® Green	26.96
RPII	TransStart® TIP Green	17.93	
	TransStart® TOP Green	20.72	
	Power SYBR® Green	36.00	

Specifications

- EvaGreen = enhanced sensitivity
- Triple HotStart blocking
- Fast Protocol Compatible
- Low /High ROX included
- 2x concentration

AQ141-01	1 ml	x2	83\$
AQ141-02	5 ml	x2	328\$
AQ141-03	15 ml	x2	837\$
AQ141-04	25 ml	x2	1174\$



**Performance
 Garantie**

À l'achat de TIP Green qPCR Supermix,
 obtenez **50% de rabais** sur l'ensemble de transcription inverse
TransScript tout-en-un! Ou obtenez 20% de rabais sans autre achat.

-50%

Reverse Transcriptase for qPCR
 TransScript All-in-One for qPCR
 (+ gDNA remover)

AT341-01	50 rx (20 ul)	Reg. 267\$ (-20%)
AT341-02	100 rx (20 ul)	Reg. 481\$ (-20%)

- gDNA remover inclus.
- Inactivation en 5 secondes!
- Inclus : RNase inhibitor et oligos
- RNase H déficient
- cDNA prêt pour le qPCR en 15 min.
- Ajoutez seulement l'ARN, H2O, gDNAr (opt.)

Promos Printemps 2018

Fin : 30 juin 2018



**NOS MEILLEURS
VENDEURS! !**

15 - 30%

DNA Polymerases

EasyTaq 10 000u &+ (-dNTPs / +dNTPs)

AP111-01 / -11	Reg. 50\$ / 58\$	-15%
AP111-02 / -12	Reg. 244\$ / 305\$	-15%
AP111-03 / -13	Reg. 693\$ / 975\$	-15%
AP111-04 (50 000u)	Reg. 1993\$	-30%

HiFi DNA Polymerases -20%

FastPfu	AP221	Reg. 143\$ - 1234\$
FastPfu FLY	AP231	Reg. 166\$ - 1375\$
KD Plus	AP301	Reg. 121\$ - 977\$

DNA Ladders

ex: 1kb, 1kb Plus, 100bp, Trans15K
à partir de 36\$/ 500ul (prix rég.)

2 pour 1

Comp. Cells + TOPO

Trans5a (DH5a), Trans10 (Top10), Trans2-Blue (XL10-Gold), TransStbl3 (Stbl3), Trans1-T1, Trans110 (dam-/dem-), pEASY Blunt, T/A cloning
à partir de 99\$ pour 5 x 100 ul

2 pour 1

Enzymes

T4 DNA Ligase, DMT (DpnI), DNase I, AvrII, BamHI, EagI, EcoRI, EcoRV, HindIII, KpnI, NcoI, NheI, NotI, PstI, SacI, SalI, ScaI, SpeI, SphI, XbaI, XhoI, XmaI

- 25 %

à partir de 35\$

Mouse Genotyping PCR Kit



- 13 min setup.
- **NO columns & NO centrifugation.**
- New **Inhibitor-resistant** Hotstart DNA Polymerase Mastermix with cyan dye.
- Multiplex PCR with up to 5 pairs of primers.

Cat.No	Format	Price (CAD)
AD501-01	100 rx	168 \$
AD501-02	500 rx	677 \$

**110\$
444\$**

MagicPure Magnetic Rack

- 16 magnets for 1.5/2.0 ml tubes





Cat.No	Format	Price (CAD)
GQ201-01	1 rack	150 \$

110\$

TransNGS
Next-Generation Sequencing
for Illumina®

- 15 %

Tn5 DNA Library Prep Kit (for 50 ng DNA) Single Index Primers
DNA Library Prep Kit for Illumina Dual Index Primers
Library Amplification SuperMix rRNA Depletion Kit (Human/Mouse/Rat)
Library Quantification DNA Standards DNA and RNA magnetic beads


Cashback x2 Promotion

 Ends June 30th, 2018

Favorgen is located in the national biotechnology park in Taiwan. Its factory is ISO-certified and possess an automated system for production of high-quality molecular biology products and clinical chemical reagents.

DNA/RNA Purification & Extraction Kits

Plasmid DNA, Gel & PCR

Genomic DNA

RNA & miRNA

FAVORGEN Kit	Cat No.	Format	Price (CAD)	CashBack x2	FAVORGEN Kits	Cat No.	Format	Price (CAD)	CashBack x2	FAVORGEN Kits	Cat No.	Format	Price (CAD)	CashBack x2
FavorPrep Plasmid DNA Extraction MiniPrep Kit (< 3ml)	FAPDE-300	300 prep	144 \$	28 \$	Blood gDNA Extraction MiniPrep Kit (+ Proteinase K)	FABGK-001-2	300 prep	300 \$	60 \$	Blood/Cultured Cell Total RNA MiniPrep Kit	FABRK-001-2	300 rx	882 \$	176 \$
FavorPrep Plasmid DNA Extraction MidiPrep Kit	FAPDE-002	25 prep	232 \$	46 \$	Blood gDNA Extraction MAXIPrep Kit (+ Proteinase K)	FABGK-003	10 prep	260 \$	52 \$	Blood/Cultured Cell Total RNA MAXIPrep Kit	FABRK-003	10 rx	294 \$	58 \$
FavorPrep Plasmid DNA Extraction MidiPrep Kit (Endo-free)	FAPDE-002-EF	25 prep	249 \$	48 \$	Tissue gDNA Extraction MiniPrep Kit (+ Proteinase K)	FATGK-001-2	300 prep	592 \$	118 \$	Tissue Total RNA MiniPrep Kit	FATRK-001-2	300 rx	882 \$	176 \$
FavorFilter Plasmid DNA Extraction MidiPrep Kit	FAFTE-002-1	50 prep	495 \$	98 \$	Plant gDNA Extraction MiniPrep Kit	FAPGK-001-2	200 prep	420 \$	84 \$	Tissue Total RNA MAXIPrep Kit	FATRK-003	10 rx	294 \$	58 \$
FavorPrep Plasmid Extraction MAXIPrep PLUS Kit (EndoToxin-Free)	FAPMX-010-EF	10 prep	240 \$	48 \$	Stool DNA Isolation MiniPrep Kit (+Proteinase K)	FASTIO-01-1	100 prep	323 \$	64 \$	Plant Total RNA MiniPrep Kit	FAPRK-001-2	300 rx	735 \$	100 \$
FavorPrep Plasmid DNA Extraction MAXIPrep Kit	FAPDE-003-1	20 prep	375 \$	74 \$	Soil DNA Isolation MiniPrep Kit	FASOI-001-1	100 prep	323 \$	64 \$	Plant Total RNA MiniPrep Kit (Woody Plant)	FAPRK-003-1	100 rx	297 \$	58 \$
FavorFilter Plasmid DNA Extraction MAXIPrep Kit	FAFTE-001	10 prep	262 \$	52 \$	Fungi/Yeast gDNA Extraction MiniPrep Kit	FAFYG-001-1	100 prep	588 \$	116 \$	Viral Nucleic Acid Extraction Kit I	FAVNK-001-2	300 rx	718 \$	142 \$
FavorFilter Plasmid DNA Extraction MAXIPrep Kit (EndoToxin-Free)	FAFTE-001-EF	10 prep	280 \$	56 \$	Food DNA Extraction Kit I	FDK-1050	50 prep	420 \$	82 \$	96-well Total RNA Kit	FATRE-96004	4 plates	890 \$	178 \$
FavorPrep GEL/PCR Purification Kit	FAGCK-001-1	300 prep	164 \$	32 \$	Milk Bacterial DNA Extraction Kit	FAMBD-001	50 prep	420 \$	82 \$	96-well Viral DNA/RNA Extraction Kit	FAVRE-96004	10 plates	890 \$	178 \$
FavorPrep MicroElute GEL/PCR Purification Kit	FAEPK-001-1	200 prep	77 \$	14 \$	Plant Genomic DNA Extraction MiniPrep Kit	FAPGK-001-2	200 prep	420 \$	82 \$	miRNA Isolation Kit	FAMIK-001	100 prep	560 \$	112 \$
FavorPrep 96-well Plasmid DNA Extraction Kit	FAPWE-96004	4 plates	445 \$	88 \$	FPPE Tissue DNA Extraction MicroElute Kit	FAFFM-050	50 prep	170 \$	34 \$	Tri-RNA Reagent	FATRR-003	450 ml	658 \$	130 \$
FavorPrep 96-well GEL/PCR Purification Kit	FAPKE-96004	4 plates	445 \$	88 \$	96-well Genomic DNA Extraction Kit	FADWE-96004	4 plates	557 \$	110 \$	AfterTri-RNA Purification Kit	FAATR-001-1	200 rx	385 \$	76 \$
FavorPrep Plasmid DNA Extraction Mini Kit (Automatable - QIAcube)	FAPDE-AUTO-024	24 prep	115 \$	22 \$	96-well Viral Nucleic Acid Extraction Kit	FAVRE-96004	4 plates	890 \$	178 \$					

* CashBack rewards are applicable on subsequent purchases within the next 3 months after purchase.
See www.civicbio.com/promotions/cashback-rewards/ for details.

More Mini/Midi/Maxi/96-well kits and sizes are available.

FavorFilter

FavorFilter Plasmid DNA Extraction **MAXIPrep** Kit
10 preps / (Endo-free option)

FAFTE-001 / -001-EF
262\$ / 280\$

FavorFilter Plasmid DNA Extraction **MidiPrep** Kit
25 / 50 preps

FAFTE-002 / -002-1
325\$ / 495\$

**Faites-vous la faveur d'utiliser (ENFIN!)
des seringues à lysat très efficaces!**

Customer Review - Québec, October 25, 2017

"Starting from a 100ml culture, we obtained

3 mg of pure DNA using this kit, while we used to get

around 500µg with the kit of competitor Q Filter kits, for the same prep."

222 \$

PMINIGP300
Valid until 2018/08/31

300 preps
Plasmid DNA **MiniPrep**
FAPDE-300 (reg. 144\$)



300 preps
GEL/PCR
Purification
FAGCK-001-1 (reg. 164\$)

Clone miniprep	OD260	OD260	OD280	Dilution factor	Conc. (ng/ul)	OD260/280 ratio
MCB-411	150	355	183	20	355	1,94
MCB-421	179	430	224	20	430	1,92
MCB-422	181	438	229	20	438	1,91

Livraison Gratuite
commandes > 100\$
25\$ si < 100\$

Récompenses Cashback

Produits Illimités

Services PCR, Clonage, Mutagénèse



Custom Experimental Design
We design everything for you.
Tailored to your needs.
No need to worry about feasibility



Plasmid Maps
Complete and detailed maps
provided in Snapgene format,
Genbank or Fasta.



Sequencing Results
DNA sequencing is performed at
the Plateforme de Séquençage
du CHUL in Quebec City.



Glycerol Stocks
In contrast to other service providers,
we provide the bacteria glycerol stocks.
Whenever possible, plasmid DNA is maintained
in T1-phage and T5-phage resistant competent cells.



Fair Pricing
Discounts and/or free complementary
products are offered for large projects or
long-term projects.

www.civicbio.com/services/pcr-services/pcr-service-pricing pour plus de détails.

Les prix varient selon l'ampleur du projet, la taille des inserts, la composition en GC et la nature du plasmide ou des bactéries compétentes désirés. De 150 à 900\$ par clone en général, sauf exceptions.

Micropipette ou Rack Magnétique **GRATUIT!**

À l'achat d'au moins 1000\$ en services, obtenez gratuitement une Pipette P series ou un rack Magnétique MagicPure GQ201.

MicroPipette autoclavable

P series Variable Volume reg. 120\$

Model	Volume	Accuracy	Repeatability
P2	0.1-2 µl	±12.0-1.5%	≤6.00-0.70%
P10	1-10 µl	±2.5-1.0%	≤1.50-0.40%
P20	2-20 µl	±2.5-1.0%	≤1.50-0.30%
P100	10-100 µl	±1.8-0.8%	≤1.50-0.15%
P200	20-200 µl	±1.8-0.8%	≤0.50-0.15%
P1000	100-1000 µl	±1.5-0.8%	≤0.30-0.15%
P5000	1000-5000 µl	±1.2-0.6%	≤0.60-0.30%
P10ml	1-10 ml	±2.5-0.6%	≤0.30-0.15%

Features

- Perfect weight** : not too light or too heavy.
- Ultra smooth piston**: comfortable for your hands during repeated utilization.
- Convenient Single-handed adjustment.**
- Autoclavable** : all parts at 121°C.
- Compatible** with standard tips.

